


Anti-Pias1 Rabbit Monoclonal Antibody For Wb And Elisa - 075925

Item Number: CM0009810



Introduction

High-quality anti-PIAS1 rabbit monoclonal antibody validated for WB and ELISA. Reacts with human and mouse PIAS1, a 72kDa E3 SUMO ligase involved in transcriptional regulation and DNA repair. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	PIAS1 Rabbit mAb
Remarks/Alias	GBP; ZMIZ3; DDXBP1; GU/RH-II; PIAS1
Species	Human
GeneID (Human)	8554
GeneID	8554
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 3-101 of human PIAS1 (075925)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075925
Protein Weight	72kDa
Shipping	Ice bag