

Anti-Pias2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075928

Item Number: CM0024953



Introduction

Rabbit polyclonal antibody against human PIAS2 (E3 SUMO-protein ligase PIAS2), validated in WB, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for studying SUMOylation, transcriptional regulation, and steroid hormone signaling.

[Learn More](#)

Parameter	Value
Product Name	PIAS2 Rabbit pAb
Remarks/Alias	DIP; MIZ1; SIZ2; ARIP3; PIASX; ZMIZ4; PIAS2
Species	Human
Gene ID	9063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-621 of human PIAS2 (NP_004662.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075928
Protein Weight	68kDa
Shipping	Ice bag