

Anti-Pias3 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q9Y6X2

Item Number: CM0014774



Introduction

Rabbit polyclonal antibody targeting human PIAS3, an E3 SUMO ligase involved in STAT signaling and sumoylation. Suitable for IF/ICC and ELISA applications, validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	PIAS3 Rabbit pAb
Remarks/Alias	ZMIZ5; PIAS3
Species	Human
Gene ID	10401
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 399-628 of human PIAS3 (NP_006090.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y6X2
Protein Weight	68kDa
Shipping	Ice bag