

Anti-Pias4 Polyclonal Antibody For Wb, Elisa - Q8N2W9

Item Number: CM0027794

Introduction

Rabbit polyclonal antibody against human PIAS4, validated for WB and ELISA with cross-reactivity to mouse. Suitable for research on SUMOylation, transcriptional regulation, and DNA repair pathways.

[Learn More](#)

Parameter	Value
Product Name	PIAS4 Rabbit pAb
Remarks/Alias	PIASY; Piasg; ZMIZ6; PIAS-gamma; PIAS4
Species	Human
GeneID	51588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-260 of human PIAS4 (NP_056981.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N2W9
Protein Weight	57kDa
Shipping	Ice bag