

Anti-Pibf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wxw3

Item Number: CM0026500

Introduction

Rabbit polyclonal antibody against human PIBF1 (PIBF/CEP90/JBTS33) for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studies on ciliogenesis, progesterone signaling, and reproductive immunology.

[Learn More](#)

Parameter	Value
Product Name	PIBF1 Rabbit pAb
Remarks/Alias	PIBF; CEP90; JBTS33; C13orf24; PIBF1
Species	Human
GeneID (Human)	10464
GeneID	10464
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PIBF1 (Q8WXW3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXW3
Protein Weight	90kDa
Shipping	Ice bag