


Anti-Pick1 Polyclonal Antibody For Wb, Elisa - Q9Nrd5

Item Number: CM0015434



Introduction

High specificity anti-PICK1 polyclonal antibody validated for Western blot and ELISA applications, targeting human PICK1, a protein involved in synaptic plasticity and receptor trafficking. Ideal for neuroscience research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PICK1 Rabbit pAb
Remarks/Alias	PICK; PRKCABP; PICK1
Species	Human
GeneID (Human)	9463
GeneID	9463
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PICK1 (NP_036539.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NRD5
Protein Weight	47kDa
Shipping	Ice bag