


Anti-Pick1 Rabbit Polyclonal Antibody For Ihc-P, Elisa - Q9Nrd5

Item Number: CM0017922



Introduction

Rabbit polyclonal antibody against PICK1, validated for IHC-P and ELISA. Suitable for investigating PRKCA-binding protein's role in synaptic plasticity and receptor trafficking in human samples.

[Learn More](#)

Parameter	Value
Product Name	PICK1 Rabbit pAb
Remarks/Alias	PICK; PRKCABP; PICK1
Species	Human
Gene ID (Human)	9463
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human PICK1 (NP_036539.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9NRD5
Protein Weight	47kDa
Shipping	Ice bag