

Anti-Pid1 Polyclonal Antibody For Wb, Ip, Elisa - Q7Z2X4

Item Number: CM0018412

Introduction

Anti-PID1 rabbit polyclonal antibody validated for WB, IP, and ELISA. Detects PID1 in human, mouse, and rat samples. Cytoplasmic protein involved in preadipocyte proliferation, suitable for obesity and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	PID1 Rabbit pAb
Remarks/Alias	PCL11; NYGGF4; P-CLI1; HMFN2073; PID1
Species	Human
Gene ID (Human)	55022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-248 of human PID1 (NP_060403.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z2X4
Protein Weight	28kDa
Shipping	Ice bag