


Anti-Pidd Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Hb75

Item Number: CM0025490



Introduction

PIDD rabbit polyclonal antibody, validated for WB, IHC-P, IF/ICC, ELISA applications; reacts with human, mouse, rat. Targets p53-induced death domain protein 1, key in DNA damage response and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	PIDD Rabbit pAb
Remarks/Alias	LRDD; PIDD; MRT75
Species	Human
GeneID (Human)	55367
GeneID	55367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human PIDD (NP_665893.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB75
Protein Weight	100kDa
Shipping	Ice bag