


Anti-Piezo1 (Fam38A) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92508

Item Number: CM0026309



Introduction

Rabbit polyclonal antibody against PIEZO1 (FAM38A), a mechanosensitive ion channel component. Validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Suitable for touch/pain and vascular research.

[Learn More](#)

Parameter	Value
Product Name	FAM38A/PIEZO1 Rabbit pAb
Remarks/Alias	ER; DHS; Mib; LMPH3; FAM38A; LMPHM6; FAM38A/PIEZO1
Species	Human
GeneID (Human)	9780
GeneID	9780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2230-2420 of human FAM38A/PIEZO1 (NP_001136336.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92508
Protein Weight	287kDa
Shipping	Ice bag