


Anti-Piezo1/Fam38A Polyclonal Antibody For Wb, Elisa - Q92508

Item Number: CM0029044



Introduction

Rabbit polyclonal antibody against human PIEZO1 mechanosensitive ion channel, validated for Western blot and ELISA. Suitable for mechanotransduction, touch sensation, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	FAM38A/PIEZO1 Rabbit pAb
Remarks/Alias	ER; DHS; Mib; LMPH3; FAM38A; LMPHM6
Species	Human
GeneID (Human)	9780
GeneID	9780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1811-1952 of human FAM38A/PIEZO1 (NP_001136336.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92508
Protein Weight	287kDa
Shipping	Ice bag