

Anti-Pif4 Rabbit Polyclonal Antibody For Wb And Elisa - Q8W2F3

Item Number: CM0019368

Introduction

High-quality rabbit polyclonal antibody against Arabidopsis thaliana transcription factor PIF4, validated for Western blot and ELISA. Ideal for phytochrome signaling, light and temperature stress response, and plant thermomorphogenesis research.

[Learn More](#)

Parameter	Value
Product Name	PIF4 Rabbit pAb
Remarks/Alias	AtPIF4; MFL8.13; MFL8_13; phytochrome interacting factor 4; SRL2; PIF4
Species	Arabidopsis thaliana
GeneID (human)	818903
GeneID	818903
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of arabidopsis thaliana PIF4 (NP_565991.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q8W2F3
Protein Weight	48kDa
Shipping	Ice bag