

Anti-Pig Cd27 Monoclonal Antibody For Fc - L7Qhb1

Item Number: CM0011819

Introduction

-conjugated rabbit anti-pig CD27 monoclonal antibody for flow cytometry. Recognizes pig CD27 antigen, a costimulatory immune-checkpoint receptor. Ideal for T-cell, NK-cell, and B-cell research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD27 mAb
Species	Pig
GeneID (Human)	939
GeneID	100520023
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-180 of pig CD27 (NP_001265692.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
SWISS	L7QHB1
Protein Weight	27kDa
Shipping	Ice bag