

Anti-Pigm Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H3S5

Item Number: CM0024267

Introduction

High-quality rabbit polyclonal antibody targeting human PIGM, appropriate for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects a band of approximately 49 kDa. Suitable for research on GPI-anchor biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	PIGM Rabbit pAb
Remarks/Alias	GPI-MT-I; PIGM
Species	Human
GeneID (Human)	93183
GeneID	93183
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-96 of human PIGM (NP_660150.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H3S5
Protein Weight	49kDa
Shipping	Ice bag