

Anti-Pigo Polyclonal Antibody For Wb And Elisa - Q8Teq8

Item Number: CM0019959

Introduction

Rabbit polyclonal antibody against human PIGO, validated for WB and ELISA, cross-reacts with mouse. Suitable for GPI-anchor biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PIGO Rabbit pAb
Remarks/Alias	HPMRS2; PIGO
Species	Human
GeneID (Human)	84720
GeneID	84720
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 950-1050 of human PIGO (NP_116023.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TEQ8
Protein Weight	119kDa
Shipping	Ice bag