

Anti-Pigp Rabbit Polyclonal Antibody For Wb And Elisa - P57054

Item Number: CM0012902

Introduction

Rabbit polyclonal antibody targeting human PIGP, a key enzyme in GPI-anchor biosynthesis. Validated for WB and ELISA, with mouse cross-reactivity. Suitable for studying GPI biosynthesis and neurodevelopmental disorders.

[Learn More](#)

Parameter	Value
Product Name	PIGP Rabbit pAb
Remarks/Alias	DCRC; DSRC; DEE55; DSCR5; PIG-P; DCRC-S; EIEE55; PIGP
Species	Human
Gene ID (Human)	51227
Gene ID	51227
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-158 of human PIGP (NP_710148.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P57054
Protein Weight	18kDa
Shipping	Ice bag