

Anti-Pigq Rabbit Polyclonal Antibody For Wb, Elisa - Q9Brb3

Item Number: CM0020925

Introduction

PIGQ Rabbit Polyclonal Antibody validated for WB, ELISA, reacts with Human, Mouse. Targets PIGQ, a subunit of the GPI-GnT complex essential for GPI-anchor biosynthesis. Ideal for GPI pathway research.

[Learn More](#)

Parameter	Description
Product Name	PIGQ Rabbit pAb
Remarks/Alias	GPI1; DEE77; EIEE77; MCAHS4; GPIBD19; c407A10.1; PIGQ
Species	Human
GeneID (Human)	9091
GeneID	9091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PIGQ (NP_683721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9BRB3
Protein Weight	84kDa
Shipping	Ice bag