


Anti-Pigr Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01833

Item Number: CM0029065



Introduction

Rabbit polyclonal antibody against human PIGR, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for mucosal immunity and IgA/IgM transcytosis research.

[Learn More](#)

Parameter	Specification
Product Name	PIGR Rabbit pAb
Remarks/Alias	PIGR
Species	Human
Gene ID	5284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-600 of human PIGR (NP_002635.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01833
Protein Weight	83kDa
Shipping	Ice bag