

Anti-Pigu Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H490

Item Number: CM0019001



Introduction

PIGU Rabbit Polyclonal Antibody (Q9H490) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat PIGU, a GPI-anchor transamidase component in ER. Ideal for GPI biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	PIGU Rabbit pAb
Remarks/Alias	GAB1; PIG-U; NEDB55; CDC91L1; GPIBD21; PIGU
Species	Human
Gene ID (Human)	128869
Gene ID	128869
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PIGU (NP_536724.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H490
Protein Weight	50kDa
Shipping	Ice bag