

Anti-Pigx Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tbf5

Item Number: CM0020836

Introduction

Rabbit polyclonal antibody against human PIGX (PIG-X), validated for WB and ELISA with rat cross-reactivity. Suitable for GPI-anchor biosynthesis studies, stabilizing subunit of mannosyltransferase I complex.

[Learn More](#)

Parameter	Value
Product Name	PIGX Rabbit pAb
Remarks/Alias	PIG-X; PIGX
Species	Human
GeneID (Human)	54965
GeneID	54965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-220 of human PIGX (NP_060331.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8TBF5
Protein Weight	29kDa
Shipping	Ice bag