

Anti-Pik3Ca Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42336

Item Number: CM0013294



Introduction

PIK3CA Rabbit pAb is a polyclonal antibody targeting human, mouse, rat PIK3CA (PI3K alpha catalytic subunit). Validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for cancer and PI3K/AKT signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	PIK3CA Rabbit pAb
Remarks/Alias	MCM; CCM4; CWS5; MCAP; PI3K; CLAPO; CLOVE; MCMTC; PI3K-alpha; p110-alpha; PIK3CA
Species	Human
Gene ID (Human)	5290
Gene ID	5290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-600 of human PIK3CA (NP_006209.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42336
Protein Weight	124kDa
Shipping	Ice bag