

Anti-Pik3Cg Polyclonal Antibody For Ihc-P, Elisa - P48736

Item Number: CM0026139



Introduction

Rabbit polyclonal antibody targeting human PIK3CG (PI3Kγ catalytic subunit), validated for IHC-P and ELISA. Ideal for research on immune cell signaling, inflammation, and cardiac function. Reliable tool for PI3K pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PIK3CG Rabbit pAb
Remarks/Alias	PI3K; PIK3; IMD97; PI3CG; PI3Kgamma; p110gamma; p120-PI3K; PIK3CG
Species	Human
Gene ID (Human)	5294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-330 of human PIK3CG (NP_002640.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P48736
Protein Weight	126kDa
Shipping	Ice bag