

Anti-Pik3Cg Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P48736

Item Number: CM0029775



Introduction

Rabbit polyclonal antibody targeting PIK3CG (PI3K gamma), suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Recombinant protein immunogen. Ideal for studying PI3K signaling and immune cell function.

[Learn More](#)

Parameter	Value
Product Name	PIK3CG Rabbit pAb
Remarks/Alias	PI3K; PIK3; IMD97; PI3CG; PI3Kgamma; p110gamma; p120-PI3K; PIK3CG
Species	Human
Gene ID (Human)	5294
Gene ID	5294
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48736
Protein Weight	126kDa
Shipping	Ice bag