

Anti-Pik3Ip1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Fe7

Item Number: CM0022501

Introduction

High-specificity rabbit polyclonal antibody targeting human PIK3IP1 (HGFL), a negative regulator of hepatic PI3K signaling. Validated for WB and ELISA; ideal for PI3K pathway studies and assay development.

[Learn More](#)

Parameter	Value
Product Name	PIK3IP1 Rabbit pAb
Remarks/Alias	HGFL; Trip; hHGFL(S); PIK3IP1
Species	Human
GeneID (Human)	113791
GeneID	113791
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human PIK3IP1 (NP_443112.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96FE7
Protein Weight	28kDa
Shipping	Ice bag