

Anti-Pik3R3/P55Pik Monoclonal Antibody For Wb - Q92569

Item Number: CM0012253

Introduction

Rabbit monoclonal antibody targeting PIK3R3/p55PIK (p55-GAMMA). Validated for WB, ELISA; cross-reacts with human, mouse, rat. Ideal for insulin signaling and PI3K pathway research.

[Learn More](#)

Parameter	Value
Product Name	PIK3R3/p55PIK Rabbit mAb
Remarks/Alias	p55; p55PIK; p55-GAMMA; PIK3R3/p55PIK
Species	Human
GeneID (human)	8503
GeneID	8503
Immunogen	Synthetic peptide within amino acids 362-461 of human PIK3R3/p55PIK (Q92569)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92569
Protein Weight	54kDa
Shipping	Ice bag