

Anti-Pik3R3/P55Pik Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92569

Item Number: CM0022102

Introduction

Rabbit polyclonal antibody against human PIK3R3/p55PIK, validated for WB, IF/ICC, ELISA. Detects endogenous p55PIK in human, mouse, and rat. Ideal for PI3K signaling and insulin pathway research.

[Learn More](#)

Parameter	Value
Product Name	PIK3R3/p55PIK Rabbit pAb
Remarks/Alias	p55; p55PIK; p55-GAMMA; PIK3R3/p55PIK
Species	Human
GeneID (Human)	8503
GeneID	8503
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92569
Protein Weight	54kDa
Shipping	Ice bag