

Anti-Pik3R3/P55Pik Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92569

Item Number: CM0018875

Introduction

Rabbit polyclonal antibody targeting human PIK3R3/p55PIK. Validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human and mouse PIK3R3. Ideal for studying PI3K signaling in insulin and cancer pathways.

[Learn More](#)

Parameter	Specification
Product Name	PIK3R3/p55PIK Rabbit pAb
Remarks/Alias	p55; p55PIK; p55-GAMMA; PIK3R3/p55PIK
Species	Human
Gene ID (Human)	8503
Gene ID	8503
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 300-388 of human PIK3R3/p55PIK
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92569
Protein Weight	54kDa
Shipping	Ice bag