

Anti-Pik3R4/Vps15 Monoclonal Antibody For Wb, If/Icc, Elisa - Q99570

Item Number: CM0011517



Introduction

Rabbit monoclonal antibody specific to human PIK3R4/VPS15 (p150). Suitable for WB, IF/ICC, ELISA with cross-reactivity to mouse. Useful for autophagy and endocytosis research.

[Learn More](#)

Parameter	Value
Product Name	PIK3R4/VPS15 Rabbit mAb
Remarks/Alias	p150; VPS15; PIK3R4/VPS15
Species	Human
GeneID (Human)	30849
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1259-1358 of human PIK3R4/VPS15 (Q99570).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99570
Protein Weight	153kDa
Shipping	Ice bag