

Anti-Pik3R4/Vps15 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99570

Item Number: CM0023009



Introduction

Rabbit polyclonal antibody targeting human PIK3R4/VPS15 (p150, PI3-kinase p150 subunit). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for autophagy, endosomal trafficking, and kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	PIK3R4/VPS15 Rabbit pAb
Remarks/Alias	p150; VPS15; PIK3R4/VPS15
Species	Human
Gene ID	30849
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human PIK3R4/VPS15 (NP_055417.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99570
Protein Weight	153kDa
Shipping	Ice bag