

Anti-Pikfyve Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2I7

Item Number: CM0025162



Introduction

Rabbit polyclonal antibody against human PIKfyve, the sole PI(3,5)P2 synthesizing enzyme. Validated for WB and ELISA in human, mouse, and rat. Ideal for endosomal dynamics, lysosomal trafficking, and SARS-CoV-2 entry research.

[Learn More](#)

Parameter	Specification
Product Name	PIKfyve Rabbit pAb
Remarks/Alias	CFD; FAB1; HEL37; PIP5K; PIP5K3; ZFYVE29; PIKfyve
Species	Human
Gene ID	200576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-545 of human PIKfyve (NP_055855.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2I7
Protein Weight	237kDa
Shipping	Ice bag