

Anti-Pilra Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukj1

Item Number: CM0022506



Introduction

Rabbit polyclonal antibody targeting human PILRA, validated for Western blotting (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for immune regulation, herpes simplex virus entry, and paired receptor research.

[Learn More](#)

Parameter	Value
Product Name	PILRA Rabbit pAb
Remarks/Alias	FDF03; PILRA
Species	Human
GeneID (Human)	29992
GeneID	29992
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-180 of human PILRA (NP_840056.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UKJ1
Protein Weight	34kDa
Shipping	Ice bag