

Anti-Pim1 Polyclonal Antibody For Wb And Elisa - P11309

Item Number: CM0017805



Introduction

Rabbit polyclonal antibody against human PIM1, validated for Western blot and ELISA applications. Cross-reacts with rat. PIM1 is a proto-oncogenic kinase critical for cell survival and proliferation. Ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	PIM1 Rabbit pAb
Remarks/Alias	PIM; PIM1
Species	Human
GeneID (Human)	5292
GeneID	5292
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human PIM1 (NP_002639.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P11309
Protein Weight	36kDa
Shipping	Ice bag