

Anti-Pim1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11309

Item Number: CM0002886



Introduction

PIM1 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat PIM1 kinase, a proto-oncogene involved in cell survival and proliferation, ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	PIM1 Rabbit mAb
Remarks/Alias	PIM; PIM1
Species	Human
Gene ID (Human)	5292
Gene ID	5292
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PIM1 (P11309).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11309
Protein Weight	36kDa
Shipping	Ice bag