

Anti-Pim1 Rabbit Polyclonal Antibody For Wb & Elisa - P11309

Item Number: CM0020573



Introduction

Rabbit polyclonal antibody targeting human PIM1 kinase, suitable for Western blot and ELISA. Detects endogenous PIM1 in human, mouse, and rat samples. Essential for cancer and cell signaling research.

[Learn More](#)

Specification	Value
Product Name	PIM1 Rabbit pAb
Remarks/Alias	PIM1; PIM; serine/threonine-protein kinase pim-1
Species	Human
Gene ID	5292
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 213-313 of human PIM1 (NP_002639.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11309
Protein Weight	35kDa/45kDa
Shipping	Ice bag