

Anti-Pim1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11309

Item Number: CM0016734



Introduction

Targets human PIM1 (36 kDa), a proto-oncogene serine/threonine kinase. Recombinant immunogen (aa 1-313). Validated for WB, IHC-P, ELISA. Cross-reactivity: human, mouse, rat. Supports oncology and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PIM1 Rabbit pAb
Remarks/Alias	PIM; PIM1
Species	Human
Gene ID (Human)	5292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-313 of human PIM1 (NP_001230115.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11309
Protein Weight	36kDa
Shipping	Ice bag