


## Anti-Pim2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9P1W9

Item Number: CM0012074



### Introduction

PIM2 Rabbit Monoclonal Antibody is validated for WB and ELISA, recognizing human, mouse, and rat PIM2. Ideal for studying proto-oncogenic signaling pathways, cell survival, and proliferation in cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PIM2 Rabbit mAb                                                                                                                      |
| Remarks/Alias    | PIM2; serine/threonine-protein kinase pim-2                                                                                          |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 11040                                                                                                                                |
| GeneID           | 11040                                                                                                                                |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 10-120 of human PIM2 (NP_006866.2) |
| Source           | Rabbit                                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                                |
| Application      | WB, ELISA                                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                    |
| SWISS            | Q9P1W9                                                                                                                               |
| Protein Weight   | 34kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                              |