


Anti-Pim2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P1W9

Item Number: CM0019841



Introduction

Rabbit polyclonal antibody targeting human PIM2, a serine/threonine kinase involved in cell survival and proliferation. Suitable for WB and ELISA with cross-reactivity to human and mouse; ~34kDa.

[Learn More](#)

Parameter	Value
Product Name	PIM2 Rabbit pAb
Remarks/Alias	PIM2
Species	Human
Gene ID (Human)	11040
Gene ID	11040
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human PIM2 (NP_006866.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P1W9
Protein Weight	34kDa
Shipping	Ice bag