

Anti-Pim3 Rabbit Polyclonal Antibody For Wb And Elisa - Q86V86

Item Number: CM0023431

Introduction

PIM3 rabbit polyclonal antibody validated for Western blot and ELISA. Targets human PIM3 (pim-3), a 36 kDa serine/threonine kinase proto-oncogene involved in apoptosis and tumorigenesis. SWISS Q86V86.

[Learn More](#)

Parameter	Value
Product Name	PIM3 Rabbit pAb
Remarks/Alias	pim-3; PIM3
Species	Human
Gene ID (Human)	415116
Gene ID	415116
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 247-326 of human PIM3 (NP_001001852.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86V86
Protein Weight	36kDa
Shipping	Ice bag