

Anti-Pimreg Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bsj6

Item Number: CM0021137

Introduction

PIMREG rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for investigating mitotic regulation, chromosome segregation, and cell cycle progression in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	PIMREG Rabbit pAb
Remarks/Alias	CATS; RCS1; FAM64A; PIMREG
Species	Human
GeneID (Human)	54478
GeneID	54478
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PIMREG (NP_001182157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BSJ6
Protein Weight	27kDa
Shipping	Ice bag