


Anti-Pin1 Polyclonal Antibody For Wb And Elisa - Q9C6B8

Item Number: CM0024045



Introduction

Anti-PIN1 polyclonal antibody validated for WB and ELISA in Arabidopsis thaliana. Recognizes auxin efflux carrier component 1 (PIN1, 67 kDa), supporting auxin transport and development studies.

[Learn More](#)

Parameter	Value
Product Name	PIN1 Rabbit pAb
Remarks/Alias	ARABIDOPSIS THALIANA PIN-FORMED 1; ATPIN1; F6D5.2; F6D5_2; PIN-FORMED 1; PIN1
Species	Arabidopsis thaliana
GeneID	843693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 391-480 of arabidopsis thaliana PIN1 (NP_177500.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9C6B8
Protein Weight	67kDa
Shipping	Ice bag