


Anti-Pin1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13526

Item Number: CM0013303



Introduction

Rabbit polyclonal antibody against human Pin1, validated for WB, IF/ICC, and ELISA. Detects Pin1 in human, mouse, and rat. Ideal for research on cell cycle regulation, isomerase activity, and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Pin1 Rabbit pAb
Remarks/Alias	DOD; UBL5; Pin1
Species	Human
Gene ID (Human)	5300
Gene ID	5300
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human Pin1 (NP_006212.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13526
Protein Weight	18kDa
Shipping	Ice bag