


Anti-Pin1 Rabbit Polyclonal Antibody - Q13526

Item Number: CM0023884



Introduction

Rabbit polyclonal antibody against human Pin1 (peptidyl-prolyl cis-trans isomerase NIMA-interacting 1) for WB and ELISA; cross-reacts with mouse and rat. Suitable for cell cycle studies and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Pin1 Rabbit pAb
Remarks/Alias	DOD; UBL5; Pin1
Species	Human
Gene ID (Human)	5300
Gene ID	5300
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human Pin1 (NP_006212.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13526
Protein Weight	18kDa
Shipping	Ice bag