

Anti-Pin4 (Par14/Par17) Rabbit Pab For Wb And Elisa - Q9Y237

Item Number: CM0023889



Introduction

Rabbit polyclonal antibody targeting human PIN4 (Par14/Par17). Validated for Western blotting (WB) and ELISA. Cross-reacts with human, mouse, and rat. Suitable for ribosome biogenesis and DNA-binding research.

[Learn More](#)

Parameter	Details
Product Name	PIN4 Rabbit pAb
Remarks/Alias	EPVH; PAR14; PAR17; hEPVH; hPar14; hPar17; PIN4
Species	Human
Gene ID	5303
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-156 of human PIN4 (NP_006214.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y237
Protein Weight	14kDa
Shipping	Ice bag