


Anti-Pink1 Polyclonal Antibody (Rabbit) For Wb, If/Icc, Ihc-P, Elisa - Q9Bxm7

Item Number: CM0030224



Introduction

Rabbit polyclonal antibody targeting human PINK1 kinase; validated for WB, IF/ICC, IHC-P, ELISA in human, mouse, and rat. Suitable for mitophagy, mitochondrial quality control, and Parkinson's disease research.

[Learn More](#)

Parameter	Value
Product Name	PINK1 Rabbit pAb
Remarks/Alias	BRPK; PARK6; PINK1
Species	Human
Gene ID (Human)	65018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 282-581 of human PINK1 (NP_115785.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXM7
Protein Weight	63kDa
Shipping	Ice bag