

Anti-Pip4K2A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P48426

Item Number: CM0027094



Introduction

Rabbit polyclonal antibody against human PIP4K2A, validated for WB, IHC-P, IF/ICC, ELISA applications. Cross-reacts with human, mouse, rat. Targets phosphatidylinositol 5-phosphate 4-kinase type-2 alpha, a kinase involved in phosphoinositide signaling and autophagy.

[Learn More](#)

Parameter	Value
Product Name	PIP4K2A Rabbit pAb
Remarks/Alias	PIPK; PI5P4KA; PIP5K2A; PIP5KIIA; PIP5KII-alpha; PIP4K2A
Species	Human
GeneID (human)	5305
GeneID	5305
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-350 of human PIP4K2A (NP_005019.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48426
Protein Weight	46kDa
Shipping	Ice bag