

Anti-Pip4K2B Polyclonal Antibody For Wb, Elisa - P78356

Item Number: CM0027187



Introduction

High-quality rabbit polyclonal antibody against human phosphatidylinositol 5-phosphate 4-kinase type-2 beta (PIP4K2B) validated for WB and ELISA, suitable for studying phosphoinositide metabolism and insulin signaling.

[Learn More](#)

Parameter	Value
Product Name	PIP4K2B Rabbit pAb
Remarks/Alias	PI5P4KB; PIP5K2B; PIP5KIIB; PIP5P4KB; PIP5KIIBeta; PIP4K2B
Species	Human
Gene ID (Human)	8396
Gene ID	8396
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human PIP4K2B (NP_003550.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78356
Protein Weight	47kDa
Shipping	Ice bag