


Anti-Pist/Gopc Monoclonal Antibody For Wb - Q9Hd26

Item Number: CM0003275



Introduction

Anti-PIST/GOPC monoclonal antibody for WB, IHC-P, IF/ICC, and ELISA. Validated for human, mouse, and rat samples. Targets PIST/GOPC protein (Q9HD26). Ideal for intracellular trafficking and degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PIST/GOPC Rabbit mAb
Remarks/Alias	CAL; FIG; PIST; GOPC1; dJ94G16.2; PIST/GOPC
Species	Human
Gene ID (Human)	57120
Gene ID	57120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PIST/GOPC (Q9HD26)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HD26
Protein Weight	51kDa
Shipping	Ice bag