


Anti-Pist/Gopc Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hd26

Item Number: CM0029764



Introduction

Rabbit polyclonal antibody targeting human GOPC (PIST/CAL/FIG) for WB and ELISA applications. Detects 51kDa band in human, mouse, and rat samples. Suitable for intracellular trafficking and receptor degradation studies.

[Learn More](#)

Parameter	Value
Product Name	PIST/GOPC Rabbit pAb
Remarks/Alias	CAL; FIG; PIST; GOPC1; dj94G16.2; PIST/GOPC
Species	Human
GeneID (Human)	57120
GeneID	57120
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 195-454 of human PIST/PIST/GOPC (NP_001017408.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HD26
Protein Weight	51kDa
Shipping	Ice bag