

Anti-Pit1/Pou1F1 Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P28069

Item Number: CM0012638



Introduction

Pit1/POU1F1 Rabbit monoclonal antibody validated for WB, IP, IF-P, IHC-P, ELISA. Detects human, mouse, rat POU1F1; 33kDa. Ideal for pituitary development and hormone gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	Pit1/POU1F1 Rabbit mAb
Remarks/Alias	PIT1; CPHD1; GHF-1; Pit-1; POU1F1a
Species	Human
Gene ID (Human)	5449
Gene ID	5449
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human POU1F1 (NP_000297.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P28069
Protein Weight	33kDa
Shipping	Ice bag