

Anti-Pitpna Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q00169

Item Number: CM0013938



Introduction

Affinity-purified rabbit polyclonal antibody targeting human PITPNA. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. PITPNA is a lipid transfer protein critical for PI/PC membrane transport and signaling.

[Learn More](#)

Parameter	Value
Product Name	PITPNA Rabbit pAb
Remarks/Alias	PITPN; VIB1A; HEL-S-36; PI-TPalpha; PITPNA
Species	Human
Gene ID (Human)	5306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 211-270 of human PITPNA (NP_006215.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00169
Protein Weight	32kDa
Shipping	Ice bag